

PACO09719

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

HGSNAT

**Uniprot**

Q68CP4

**Synonyms:**

heparan-alpha-glucosaminide N-acetyltransferase; HGSNAT; DKFZp686G24175; FLJ22242; FLJ32731; HGNAT; MPS3C; TMEM76 ;

**Immunogen:**

Human HGSNAT.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A