

HMBOX1 Antibody



PACO09740

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

HMBOX1

Uniprot

Q6NT76

Synonyms:

homeobox containing 1; HMBOX1; FLJ21616; HNF1LA; PBHNF ;

Immunogen:

Human HMBOX1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A