

HNF1B Antibody



PACO09756

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

HNF1B

Uniprot

P35680

Synonyms:

HNF1 homeobox B; HNF1B; FJHN; HNF1beta; HNF2; HPC11; LF-B3; LFB3; MODY5; TCF2; VHNF1 ;

Immunogen:

Human HNF1B.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A