

HNRNPA0 Antibody



PACO09758

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

HNRNPA0

Uniprot

Q13151

Synonyms:

heterogeneous nuclear ribonucleoprotein A0; HNRNPA0; HNRPA0 ;

Immunogen:

Human HNRNPA0.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A