

# HOMER3 Antibody



PACO09771

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

HOMER3

**Uniprot**

Q9NSC5

**Synonyms:**

homer homolog 3 (Drosophila); HOMER3; HOMER-3; VESL3 ;

**Immunogen:**

Human HOMER3.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A