

HOOK3 Antibody



PACO09774

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

HOOK3

Uniprot

Q86VS8

Synonyms:

hook homolog 3 (Drosophila); HOOK3; FLJ31058; HK3 ;

Immunogen:

Human HOOK3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A