

HOXC4 Antibody



PACO09778

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

HOXC4

Uniprot

P09017

Synonyms:

homeobox C4; HOXC4; HOX3; HOX3E; cp19 ;

Immunogen:

Human HOXC4.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A