

HOXD12 Antibody



PACO09779

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

HOXD12

Uniprot

P35452

Synonyms:

homeobox D12; HOXD12; HOX4H ;

Immunogen:

Human HOXD12.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A