

HPS3 Antibody



PACO09790

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

HPS3

Uniprot

Q969F9

Synonyms:

Hermansky-Pudlak syndrome 3; HPS3; DKFZp686F0413; FLJ22704; SUTAL ;

Immunogen:

Human HPS3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A