

HPS6 Antibody



PACO09792

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

HPS6

Uniprot

Q86YV9

Synonyms:

Hermansky-Pudlak syndrome 6; HPS6; FLJ22501; MGC20522; RP11-302K17.1 ;

Immunogen:

Human HPS6.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A