

# HRASLS5 Antibody



PACO09795

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

HRASLS5

**Uniprot**

Q96KN8

**Synonyms:**

Ca(2+)-independent N-acyltransferase; H-rev107-like protein 5; HRAS-like suppressor 5; HRSL5; HRASLS5; HRLP5

**Immunogen:**

Human HRASLS5.

**Storage:**

Purified Rabbit polyclonal in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

## Product Images

---

N/A

N/A