

HS1BP3 Antibody



PACO09800

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

HS1BP3

Uniprot

Q53T59

Synonyms:

HCLS1 binding protein 3; HS1BP3; ETM2; FLJ14249; HS1-BP3 ;

Immunogen:

Human HS1BP3.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A