

HSBP1 Antibody



PACO09803

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

HSBP1

Uniprot

O75506

Synonyms:

heat shock factor binding protein 1; HSBP1; DKFZp686D1664; DKFZp686O24200; NPC-A-13 ;

Immunogen:

Human HSBP1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A