

HSD17B4 Antibody



PACO09810

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:**Protein Background:****Gene ID:**

HSD17B4

Uniprot

P51659

Synonyms:

hydroxysteroid (17-beta) dehydrogenase 4; HSD17B4; DBP; MFE-2; SDR8C1 ;

Immunogen:

Human HSD17B4.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A