

HSD17B8 Antibody



PACO09812

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

HSD17B8

Uniprot

Q92506

Synonyms:

hydroxysteroid (17-beta) dehydrogenase 8; HSD17B8; D6S2245E; FABG; FABGL; H2-KE6; HKE6; KE6; RING2; SDR30C1; dJ1033B10.9 ;

Immunogen:

Human HSD17B8.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A