

HSDL1 Antibody



PACO09815

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

HSDL1

Uniprot

Q3SXM5

Synonyms:

hydroxysteroid dehydrogenase like 1; HSDL1; MGC125994; MGC125995; MGC126032; SDR12C3 ;

Immunogen:

Human HSDL1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A