

HSDL2 Antibody



PACO09816

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

HSDL2

Uniprot

Q6YN16

Synonyms:

hydroxysteroid dehydrogenase like 2; HSDL2; C9orf99; FLJ25855; MGC10940; SDR13C1 ;

Immunogen:

Human HSDL2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A