

HSF5 Antibody



PACO09820

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

HSF5

Uniprot

Q4G112

Synonyms:

heat shock transcription factor family member 5; HSF5; FLJ40311; MGC134827 ;

Immunogen:

Human HSF5.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A