

# HSPB11 Antibody



PACO09831

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:****Protein Background:****Gene ID:**

HSPB11

**Uniprot**

Q9Y547

**Synonyms:**

heat shock protein family B (small), member 11; HSPB11; C1orf41; HSPCO34; PP25 ;

**Immunogen:**

Human HSPB11.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A