

HTRA1 Antibody



PACO09840

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

HTRA1

Uniprot

Q92743

Synonyms:

HtrA serine peptidase 1; HTRA1; ARMD7; HtrA; L56; ORF480; PRSS11 ;

Immunogen:

Human HTRA1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A