

# HUS1B Antibody



PACO09844

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

HUS1B

**Uniprot**

Q8NHY5

**Synonyms:**

HUS1 checkpoint homolog b (S. pombe); HUS1B; MGC126746; MGC126748; RP11-532F6.1 ;

**Immunogen:**

Human HUS1B.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A