

HYLS1 Antibody



PACO09847

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

HYLS1

Uniprot

Q96M11

Synonyms:

hydrolethalus syndrome 1; HYLS1; FLJ32915; HLS ;

Immunogen:

Human HYLS1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A