

IARS2 Antibody



PACO09849

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

IARS2

Uniprot

Q9NSE4

Synonyms:

isoleucyl-tRNA synthetase 2, mitochondrial; IARS2; FLJ10326 ;

Immunogen:

Human IARS2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A