

ICAM1 Antibody



PACO09851

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

ICAM1

Uniprot

P05362

Synonyms:

intercellular adhesion molecule 1; Icam1; CD54; Icam-1; Ly-47; MALA-2; MGC6195 ;

Immunogen:

Human Icam1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A