

ICAM3 Antibody



PACO09853

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

ICAM3

Uniprot

P32942

Synonyms:

intercellular adhesion molecule 3; ICAM3; CD50; CDW50; ICAM-R ;

Immunogen:

Human ICAM3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A