

ICAM4 Antibody



PACO09854

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ICAM4

Uniprot

Q14773

Synonyms:

intercellular adhesion molecule 4 (Landsteiner-Wiener blood group); ICAM4; CD242; LW ;

Immunogen:

Human ICAM4.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A