

# ICK Antibody



PACO09856

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

ICK

**Uniprot**

Q9UPZ9

**Synonyms:**

intestinal cell (MAK-like) kinase; ICK; ECO; KIAA0936; LCK2; MGC46090; MRK ;

**Immunogen:**

Human ICK.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A