

IFFO1 Antibody



PACO09870

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:**Protein Background:****Gene ID:**

IFFO1

Uniprot

Q0D2I5

Synonyms:

intermediate filament family orphan 1; IFFO1; DKFZp586I2223; FLJ20703; HOM-TES-103; IFFO; MGC117359 ;

Immunogen:

Human IFFO1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A