

IFIH1 Antibody



PACO09872

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

IFIH1

Uniprot

Q9BYX4

Synonyms:

interferon induced with helicase C domain 1; IFIH1; Hlcd; IDDM19; MDA-5; MDA5; MGC133047 ;

Immunogen:

Human IFIH1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A