

IFITM1 Antibody



PACO09875

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

IFITM1

Uniprot

P13164

Synonyms:

interferon induced transmembrane protein 1 (9-27); IFITM1;9-27; CD225; IFI17; LEU13 ;

Immunogen:

Human IFITM1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A