

IFITM2 Antibody



PACO09876

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

IFITM2

Uniprot

Q01629

Synonyms:

interferon induced transmembrane protein 2 (1-8D); IFITM2;1-8D ;

Immunogen:

Human IFITM2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A