

IFT140 Antibody



PACO09885

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

IFT140

Uniprot

Q96RY7

Synonyms:

intraflagellar transport 140 homolog (Chlamydomonas); IFT140; DKFZp564L232; FLJ10306; FLJ30571; KIAA0590; WDTC2; c305C8.4; c380F5.1; gs114 ;

Immunogen:

Human IFT140.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A