

IFT20 Antibody



PACO09886

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

IFT20

Uniprot

Q8IY31

Synonyms:

intraflagellar transport 20 homolog (Chlamydomonas);

Immunogen:

Human IFT20.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A