

IFT57 Antibody



PACO09888

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

IFT57

Uniprot

Q9NWB7

Synonyms:

intraflagellar transport 57 homolog (Chlamydomonas); IFT57; ESRRBL1; FLJ10147; HIPPI; MHS4R2 ;

Immunogen:

Human IFT57.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A