

# IFT88 Antibody



PACO09890

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF, IP

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

IFT88

**Uniprot**

Q13099

**Synonyms:**

intraflagellar transport 88 homolog (Chlamydomonas); IFT88; D13S1056E; DAF19; MGC26259; RP11-172H24.2; TG737; TTC10; hTg737 ;

**Immunogen:**

Human IFT88.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A