

# IGKV1-5 Antibody



PACO09898

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

IGKV1-5

**Uniprot**

P01602

**Synonyms:**

immunoglobulin &kappa; variable 1-5; IGKV1-5; IGKV; IGKV15; L12; L12a; MGC22745; MGC32715; MGC88810; V1 ;

**Immunogen:**

Human IGKV1-5.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A