

IGSF11 Antibody



PACO09899

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

IGSF11

Uniprot

Q5DX21

Synonyms:

immunoglobulin superfamily, member 11; IGSF11; BT-IgSF; CXADRL1; Igsf13; MGC35227; VSIG3 ;

Immunogen:

Human IGSF11.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A