

IGSF21 Antibody



PACO09900

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

IGSF21

Uniprot

Q96ID5

Synonyms:

immunoglobulin superfamily, member 21; IGSF21; FLJ41177; MGC15730 ;

Immunogen:

Human IGSF21.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A