

# IGSF8 Antibody



PACO09901

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

IGSF8

**Uniprot**

Q969P0

**Synonyms:**

immunoglobulin superfamily, member 8; IGSF8; CD316; CD81P3; EWI2; PGRL ;

**Immunogen:**

Human IGSF8.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A