

IL1B Antibody



PACO09922

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

IL1B

Uniprot

P01584

Synonyms:

interleukin 1, beta; IL1B; IL-1; IL1-BETA; IL1F2 ; IL1beta

Immunogen:

Human IL1B.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A