

# IL3 Antibody



PACO09934

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

IL3

**Uniprot**

P08700

**Synonyms:**

interleukin 3 (colony-stimulating factor, multiple); IL3; IL-3; MCGF; MGC79398; MGC79399; MULTI-CSF ;

**Immunogen:**

Human IL3.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A