

IL7R Antibody



PACO09939

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

IL7R

Uniprot

P16871

Synonyms:

interleukin 7 receptor; IL7R; CD127; CDW127; IL-7R-alpha; IL7RA; ILRA ;

Immunogen:

Human IL7R.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A