

CXCR2 Antibody



PACO09940

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

CXCR2

Uniprot

P25025

Synonyms:

interleukin 8 receptor, beta; IL8RB; CD182; CDw128b; CMKAR2; CXCR2; IL8R2; IL8RA ;

Immunogen:

Human IL8RB.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A