

PACO09946

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

IMMT

Uniprot

Q16891

Synonyms:

inner membrane protein, mitochondrial (mitofilin); IMMT; DKFZp779P1653; HMP; MGC111146; P87; P87/89; P89; PIG4; PIG52 ;

Immunogen:

Human IMMT.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A