

# IMP4 Antibody



PACO09948

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:****Protein Background:****Gene ID:**

IMP4

**Uniprot**

Q96G21

**Synonyms:**

IMP4, U3 small nucleolar ribonucleoprotein, homolog (yeast); IMP4; BXDC4; MGC19606 ;

**Immunogen:**

Human IMP4.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A