

# INE1 Antibody



PACO09951

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

INE1

**Uniprot**

O15225

**Synonyms:**

inactivation escape 1 (non-protein coding); INE1; MGC97326; MGC97335; MGC97386; NCRNA00010 ;

**Immunogen:**

Human INE1.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A