

# INO80 Antibody



PACO09960

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

INO80

**Uniprot**

Q9ULG1

**Synonyms:**

INO80 homolog (S. cerevisiae); INO80; INO80A; INOC1; hINO80 ;

**Immunogen:**

Human INO80.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A