

INPP5J Antibody



PACO09964

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

INPP5J

Uniprot

Q15735

Synonyms:

inositol polyphosphate-5-phosphatase J; INPP5J; INPP5; MGC129984; PIB5PA; PIPP ;

Immunogen:

Human INPP5J.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A