

# INTS12 Antibody



PACO09971

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

INTS12

**Uniprot**

Q96CB8

**Synonyms:**

integrator complex subunit 12; INTS12; INT12; PHF22; SBB122 ;

**Immunogen:**

Human INTS12.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A