

INTS8 Antibody



PACO09976

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

INTS8

Uniprot

Q75QN2

Synonyms:

integrator complex subunit 8; INTS8; C8orf52; FLJ20530; INT8; MGC131633 ; INT8

Immunogen:

Human INTS8.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A